

- ✓ Improve Accuracy and Consistency
- ✓ Increase Convenience and Productivity

0.5% Methyl Cellulose Kit



Use this concentration for the wash cycle between injections. Filtered to ensure consistent viscosity to coat the capillary lumen to minimize EOF. Conveniently packaged in 2 x 100mL containers. Expiry dating is 12 months. Reorder No: 102505

Electrolyte Kit



0.08M Phosphoric Acid Solution H_3PO_4
0.1M Sodium Hydroxide Solution $NaOH$

These 100mL Analyte and Catholyte solutions in 0.1% MC are used to fill the electrolyte tanks on the iCE280 cartridge. The labels and bottles are color coded to improve safety. Expiry dating is 12 months. Reorder No: 102506

Chemical Test Kit



The Chemical Test Kit is designed to confirm the overall performance of the iCE280 on 3 separate occasions. Each time, performance can be checked repeatedly over a 3 day period. The kit contains anolyte, catholyte, MC and a standard sample with Hemoglobin AFSC. Expiry dating is 6 months. Reorder #: 101801F

System Suitability Kit



Use the System Suitability Kit to run a control test of the iCE280 on 10 separate occasions. Each time, the performance can be checked up to 6 times over 3 days. The kit contains a standard sample with Hemoglobin AFSC. Expiry dating is 6 months. Reorder #: 102093F



To Order
toll free 1.800.206.1027
fax 416.231.1152
sales@convergentbiosci.com

Nov 10 08